

A New Approach For Co-localization Of Proteins Using HRP And AP Enzyme Detection In Paraffin Embedded Tissues

Rachel M Gonzalez, Li Ding, Ming Xiao, Xiaoping Zhang
GBI Labs-Inc, Mukilteo, WA 98275 info@gbi-inc.com



Abstract Co-localization of two proteins is a powerful tool in immunofluorescence, however the stability of the fluorescence limits the ability to screen large sample sets. Double and triple staining with the ability to visualize co-localization of proteins using HRP and AP enzymes in paraffin embedded tissue may prove to be a useful tool for discovery and functional analysis of gene and protein expression. Here we present a novel chromogen mix by GBI Labs Inc. of Emerald and Red that yields a third color (Dark Blue/Purple) when proteins are co-localized using HRP and AP enzymes for detection. Initial tests evaluated combinations of ER, Her2, CEA, Keratin, pR, p53, PCNA, and Ki67 in 10 breast and 10 colon cancer cases to screen the percentage of co-localized proteins in these cases. In addition, a 70 case breast cancer tissue array containing duplicate samples and known diagnosis of ER, PR, and Her2 for each case were screened with ER, PR, Her2, p53, PCNA, or Ki67 and keratin antibodies. The core tumor sections showed co-localization of PCNA, Ki67 and p53 proteins. These new double and triple staining protocols in combination with the tissue array, will allow researchers to rapidly determine the association between expression of two or more proteins in the tissue and if proteins co-localize within the cell.

Staining Pattern of Proteins with AP-Red+ & Emerald Chromogen




Fig. 1 Cartoon of staining pattern of proteins expressed separately versus co-localized in the cell.

Methods Colon and breast cancer tissue specimens were obtained from Dr. Shi, Beijing, China and the 70 case breast cancer tissue array was obtained from BioChain Cat# Z7020004. Representative cores were used to generate two tissue arrays, one of breast cancer and the second of colon cancer. Screens were done as follows. Five-µm sections were cut and mounted on coated slides and dried overnight at 37C. Slides were de-waxed in xylene and hydrated using graded alcohols then rinsed in tap water. Endogenous peroxidase was blocked with 3% H2O2 for 10 minutes and washed in several changes of water. If Heat Induced Epitope Retrieval was required for primary antibody this step was done with 10mM Citrate buffer pH6.0 for 15 min at 98C with a cool-down to 45C. The double staining procedure is highly dependent on the primary antibody combination with respect to animal species. Mouse –Rabbit and Mouse –Goat primary antibody combinations were incubated on the tissue together. Mouse- Mouse and Rabbit-Rabbit primary antibody combinations were stained sequentially following protocols used in the double and triple stain kits listed below. Most primary antibody combinations were incubated for 30 minutes at room temperature unless otherwise indicated by primary antibody source. AP-Red+ chromogen requires incubation prior Emerald chromogen as detailed in the multi-staining protocol. Tissues were scored using Olympus BX40 Light microscope. Scoring was assessed on total percent positive cells and intensity of stain. Co-localization was assessed on all levels of positive cells.

Antibody Pairs	Source	Tissue	Double & Triple Staining Kts GBI Labs
Ms anti-p53 & Rb anti-ER	GBI Labs M16-08 Spring E1644	Breast Cancer	DS201C, DS202C, DS233C
Rb anti-ER & Ms anti-PCNA	Spring E1644 GBI Labs M16-01	Breast Cancer	DS203C, DS212C
Ms anti-ER & Rb anti-PR Rb anti-HER2	GBI Labs ERC Epitomics 148-1 Epitomics 420-1	Breast Cancer	TS302A, TS309A
Ms anti-p53 & Ms anti-PCNA & Rb anti-CEA	GBI Labs M16-08 GBI Labs M16-01 GBI Labs P3012	Colon cancer	TS301A, TS308A
Gt anti-EGFR Ms anti-CEA	Santa Cruz S203-G Gift from Dr. Shi, ZCB	Colon Cancer	DS207C,
Gt anti-EGFR Ms anti-CK8/18	Santa Cruz S203-G GBI Labs M308	Breast Cancer	DS207C
Ms anti-Actin Rb anti-Desm	Sigma A2547-5ml GBI Labs P4005	Colon Cancer	DS201C, DS202C, DS233C

Results




Fig. 2 Expression pattern of Rb anti-ER with Emerald and Ms anti-p53 with AP-Red+ (Panel A&D) on two breast cancers: Rb anti-ER detection with emerald chromogen (Panel B&E); and Ms anti-p53 with AP-Red+ chromogen (Panel C&F). Co-localization represented in blue/purple color. Although the stain for these two tumors show similar expression, the amount of co-localization is very different as seen in panel A with >70% of tumor cells expressing both proteins and panel D with 25% of the tumor cells expressing both proteins as indicated by blue/purple color.




Fig. 3 Two Colon Cancer double stained with Rb anti-CEA with emerald and Ms anti-CK8/18 with AP-Red+ (Panel A&D); Rb anti-CEA only detection with emerald chromogen (Panel B&E); and Ms anti-CK8/18 only with AP-Red+ chromogen. Co-localization represented in blue/purple color. Panel A shows that when tumor is negative, there is no background or overlap; only when the two proteins are co-expressed there is a color change.




Fig. 4 Panel (A) colon cancer double stained with Rb anti-Desm/Emerald and Ms anti-Smooth Muscle Actin/AP-Red+. Panel (B, C, E, & F) are breast cancer cases. Panel (B & E) are stained with Ms anti-CK8/18 with AP-Red+ and Gt anti-EGFR with emerald chromogen. Panel (C & F) are stained with Ms anti-ER with AP-Red+ and Rb anti-PR with emerald chromogen. Panel (D) colon cancer stained with Gt anti-EGFR emerald & Ms anti-CEA/AP-Red+.

Tables BioChain Breast Tumor Tissue Array (BCTA) double stain with Ms anti-PCNA & Gt anti-p53 and Ms anti-CK8/18 and Rb anti-HER2

Tissue	Ms anti-PCNA & Gt anti-p53												
	1	2	3	4	5	6	7	8	9	10			
1	NA	NA	3%	20%	10%	NA	1%	NA	NA	NA	NA	5%	NA
2	NA	NA	1%	20%	90%	NA	NA	NA	NA	10%	NA	5%	NA
3	NA	NA	NA	NA	NA	30%	10%	NA	NA	NA	30%	54%	10%
4	NA	1%	NA	NA	NA	40%	10%	10%	NA	NA	1%	64%	42%
5	NA	0.5%	NA	NA	NA	NA	50%	NA	NA	5%	NA	50%	20%
6	NA	NA	NA	NA	NA	40%	45%	1%	NA	2%	NA	50%	16%
7	NA	NA	NA	50%	NA	70%	NA	NA	NA	50%	NA	30%	NA
8	50%	NA	NA	45%	20%	36%	NA	NA	NA	10%	NA	10%	1%
9	5%	NA	NA	NA	NA	NA	60%	NA	NA	NA	NA	NA	NA
10	NA	NA	NA	NA	NA	NA	NA	50%	NA	NA	NA	NA	NA

Tissue	Rb anti-HER2 & Ms anti-CK8/18												
	1	2	3	4	5	6	7	8	9	10			
1	NA	NA	NA	5%	NA	NA	NA	20%	80%	NA	54%	NA	NA
2	NA	NA	NA	30%	NA	NA	NA	10%	64%	NA	63%	NA	NA
3	NA	NA	NA	NA	80%	70%	90%	80%	NA	20%	NA	50%	50%
4	NA	90%	NA	NA	80%	70%	90%	90%	NA	10%	NA	50%	50%
5	NA	10%	NA	NA	NA	48%	NA	NA	NA	NA	20%	NA	NA
6	10%	10%	NA	90%	NA	40%	NA	NA	NA	70%	20%	NA	NA
7	NA	NA	40%	NA	NA	NA	NA	NA	NA	NA	NA	NA	NA
8	NA	NA	40%	NA	NA	70%	NA	NA	NA	NA	NA	NA	NA
9	NA	NA	70%	NA	NA	NA	80%	NA	35%	NA	NA	NA	NA
10	25%	NA	100%	NA	NA	NA	42%	NA	25%	NA	NA	NA	NA




Fig. 5 Expression pattern of double and triple stain Rb anti-CEA with DAB, Ms anti-p53 with Emerald, and Ms anti-PCNA with AP-Red+ chromogen detection (Panel A&B) on colon cancer. Single stains are shown in panel C, D, & E.




Fig. 6 BioChain Breast tumor Tissue Array double stain with (A) Ms anti-PCNA (AP-Red+) & Gt anti-p53 (Emerald) and (B) Ms anti-CK8/18 (AP-Red+) and Rb anti-HER2 (Emerald).

Conclusion: Here we presented multiple double staining showing that the emerald chromogen in the presence of the AP-Red+ chromogen produces a blue/purple color when two proteins are co-expressed in the nucleus, cytoplasm, or cell membrane. We predict that a semi-quantitative assessment can be made based on which protein level is expressed at higher levels which results in how blue (more emerald antigen) or how purple (more red AP Red+ antigen) the stain presents. The emerald chromogen will enable the researcher to take a closer look at the co-expression of proteins in larger study sets with the use of light microscope.